


Pulse Oximeter User Manual



Product Model: YM101/YM102/YM103

Version: 2.0
Date: 2020-07-08

1 Product Introduction and Operation Guide

1.1 Front View




Figure 1 Front View of YM101/YM102/YM103

1.2 Operation Method

- A. Open the battery cover, and put the two AAA batteries into the battery compartment in correct polarities, then replace the cover.
B. Press the bottom of the equipment and open the probe, then insert one finger into the probe.
C. Press the button to turn the equipment on, and the measure interface will appear.
D. After about 8 seconds, the measurement result can be read directly from the display screen.
E. Before reading the parameters, make sure that stable numbers of the pulse oximeter interface has sustained more than 4 seconds.
F. The equipment will turned off automatically within 8 seconds when the finger left the probe.

1

1.3 Battery Installation

- A. Put the two AAA batteries into battery compartment in correct polarities (Figure 2).
B. Push the battery cover horizontally along the arrow shown as right.

WARNINGS:

- Battery polarities should be correctly installed, otherwise, damage may be caused to the equipment.
Please remove the batteries if the equipment will not use for a long time.




Figure 2 Battery Installation

1.4 Lanyard installation

- A. Pass the thinner end of the lanyard through the hanging hole;
B. Pass the thicker end of the lanyard through the thinner end and tighten the lanyard (Figure 3).

2




Figure 3 Lanyard Installation

1.5 Attention for Operation

- A. Before use check and confirm that the people or finger size were applicable.
B. Before use check and confirm that the environment should be non-combustible material, as well as to avoid high or low temperature and humidity, but also need to pay attention to the following:
a) To avoid glare and direct sunlight exposure;
b) To avoid radiation infrared or ultraviolet radiation;
c) Avoid contact with the organic solvent, mist, dust, corrosive gases.
C. The equipment should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
D. The equipment may not work normally on microcirculation barrier patients. Warm or rub the finger, or re-position the equipment could improve the measurement.
E. The ray between photo detector and light emitting diode should across patient's arteriole.
F. The patient should not use enamel or other makeup;
G. Avoid to insert a wet finger into the probe.

Figure 4 Finger Placement Diagram

1.6 Functions and Menu Operation

1.6.1 The button operation rules

- Long-press functions include entering menus, activating item's submenu, confirming setting values, and exiting item's submenu; short-press functions are polling menu items and viewing the setting values of items. It should be noted that long-press means pressing the key for about 2 seconds, and short-press means pressing the key for less than 0.5 second.

3

Notes:

- A. The user should fully insert the finger into the probe.
B. It is recommended to let the LED light shine directly on the nail (Figure 4).
C. Don't shake the finger and try to keep the patient still during the measurement.




Figure 4 Finger Placement Diagram

should be noted that long-press means pressing the key for about 2 seconds, and short-press means pressing the key for less than 0.5 second.

1.6.2 Menu Operation

Active the menu

- After the oximeter is turned on, long-press the power button to activate the menu, then short-press the button to view the setting values of each item. If the user wants to change the setting value of the item, long-press to enter the item's submenu, the parameter value starts to flash, short-press to traverse the parameter value until the parameter value required by the user is selected, long press to confirm and exit the submenu.

Item1. Setup the LED display brightness

- The first item is to setup the display brightness. Long-press the button to select a brightness level ranging from 1 to 3. The greater the value, the greater brightness of the display.



Item2. Setup the SpO2 Alarm Limits

- The second item is to setup the SpO2 alarm limits. For example: When SpO2 High limit is set to 96, an alarm will be issued when the spO2 value is higher than 96, and when


4

SpO2 low limit is set to 94, an alarm will be issued when the spO2 value is lower than 94.




Item3. Setup the PR Alarm Limits

- The third item is to setup the PR alarm limits. For example: When PR High limit is set to 130, an alarm will be issued when the PR value is higher than 130, and when PR low limit is set to 50, an alarm will be issued when the PR value is lower than 50.



Item4. Turn Alarm On/Off

- The fourth item is long-press to turn Alarm on/off.



Item5. Check the software version

- The fifth item is to view the software version.

5

2 Specifications

2.1 Classification

- Type of protection against electric shock: II (Internally powered equipment)
Degree of protection against electric shock: Type BF-Applied part
Operating mode: Shock checking
Degree of protection against hazards of explosion: IP22

2.2 Power Requirements

- Specification of battery: Two AAA (LR03)
Operating current: 25-50mA

2.3 Physical Specifications

- Width*Height*Depth: 57*30*31 mm
Weight: 28g (Bare machine)

2.4 Measurement Specifications

- SpO2 declared accuracy: 70%~100%: ±2digits
0%~60%: unspecified
SpO2 Display Range: 30%~99%
SpO2 Resolution: 1%
PR declared accuracy: 25~250 bpm: ±3digits
PR Resolution: 1 bpm

7

2.5 Environmental Specifications

- Temperature: Operating: +50~+104°F / +10~+40°C; Storage/Transportation: -4~+140°F / -20~+60°C
Humidity: Operating: 15~95%, noncondensing; Storage/Transportation: 10~95%, noncondensing
Atmosphere Pressure: Operating: 70~106kpa; Storage/Transportation: 50~107.4kpa

2.6 Display

- Display Type: 1.5" LED Display;
Display Color: YM101: Red; YM102: Green;
SpO2%, Pulse Rate, PR%, Bar Graph, Battery Indicator

Notes:

- 1) The claim for oxygen saturation accuracy should be supported by clinical studies covering the entire claimed range. The fraction of inspired oxygen (FiO2) delivered to test subjects is varied to achieve a series of targeted steady-state saturation periods over the specified SpO2 accuracy range (e.g. 70 % to 100 %), then the SpO2 accuracy is calculated by comparing SpO2 readings of the pulse oximeter to the values of SaO2 determined with a Co-Oximeter.
2) The clinical trial included 11 subjects, including 6 males and 5 females, with an age range of 18 to 46 years, the subjects skin color included dark black, medium black, light color and white.

8

3 Maintenance, Cleaning, Disinfection

3.1 Maintenance

- The equipment's design life expectancy is about 2 years, keep your equipment and accessories free of dust and dirt, and follow these rules:
A. Please clean the equipment before use according to chapter 3.2. Remove the batteries inside the battery cassette if the equipment will not be operated for a long time.
B. Replace the batteries in time when the battery voltage indicate lamps were empty.
C. It is recommended that the equipment should be kept in a dry environment with no corrosive gases and good ventilation anytime. The moisture and high-light environments will affect its lifetime and even might damage the equipment.
D. It is best to preserve the product in a place where the temperature is between -20 to 60°C and the relative humidity is less than 95%.
E. The packed equipment can be transported by ordinary conveyance. The equipment not be transported mixed with toxic, harmful, corrosive materials.

Warnings

- No modification of this equipment is allowed.

9

3.2 Cleaning

- Your equipment should be cleaned on a regular basis. If there is heavy pollution or lots of dust and sand in your place, the equipment should be cleaned more frequently. Before cleaning the equipment, consult your hospital's regulations for cleaning the equipment. Recommended cleaning agents are:
a) Mild soap (diluted);
b) Ethanol (70%);
c) Clean the exterior surface of the equipment and probe using a soft cloth dampened with the cleaner;
d) Wipe off all the cleaning solution with a dry cloth after cleaning if necessary;
e) Dry your equipment in a ventilated, cool place. To avoid damage to the equipment, follow these rules:

CAUTIONS:

- Always dilute according the manufacturer's instructions or use lowest possible concentration.
Do not immerse part of the equipment in the liquid.
Do not pour liquid onto the equipment or accessories.
Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).
If you spill liquid onto the equipment, contact us or your service personnel.

10

Certificate
PN: Finger-Pulse Oximeter
Date:

3.3 Disinfection

Clean the pulse oximeter before disinfecting it. The recommended disinfectant is ethanol 70%. Disinfection steps are the same as cleaning.

CAUTION

- Never use ETO or formaldehyde for disinfection.

3.4 Disposal

Dispose of the pulse oximeter in accordance with local environment and waste disposal laws and regulations.

4 Accessories

- One lanyard.
Two AAA batteries (Optional).
One user manual.
One certificate card.

5 Troubleshooting

Table with 3 columns: Trouble, Possible Reason, Solutions. Rows include SpO2 and PR not displayed normally, display off suddenly, SpO2 and PR not displayed stably, and device can't be turned on.

11

Table with 3 columns: Trouble, Possible Reason, Solutions. Rows include SpO2 and PR not displayed normally, display off suddenly, SpO2 and PR not displayed stably, and device can't be turned on.

12

6 Appendix A EMC

The equipment complies with the requirement of standard EN 60601-1-2:2014 "Electromagnetic Compatibility - Medical Electrical Equipment".

Table with 3 columns: Test level, Compliance level, Electromagnetic environment - guidance. Rows include emissions test, immunity test, and electrostatic discharge.

13

Table with 3 columns: Voltage fluctuations / flicker emissions, Test level, Compliance level. Includes EMC 61000-3-3 and EMC 61000-4-2 tests.

14

Table with 3 columns: Electrostatic transient / burst, Surge, Voltage dips, short interruptions and voltage variations on power supply input lines, Test level, Compliance level. Includes IEC 61000-4-4, IEC 61000-4-5, and IEC 61000-4-11 tests.

15

Table with 3 columns: Power frequency magnetic field, Power frequency magnetic fields, Conducted RF, Test level, Compliance level. Includes IEC 61000-4-8, IEC 61000-4-9, and IEC 61000-4-6 tests.

16

Table with 3 columns: Radiated RF, Test level, Compliance level. Includes IEC 61000-4-3 test.

17

Table with 3 columns: Recommended separation distance, Formulae, Test level, Compliance level. Includes IEC 61000-4-3 test.

18

Table with 3 columns: Interference may occur in the vicinity of equipment, Field strengths from fixed transmitters, Field strengths from fixed RF transmitters, Test level, Compliance level. Includes IEC 61000-4-3 test.

19

Table with 3 columns: Field strengths from fixed transmitters, Field strengths from fixed RF transmitters, Test level, Compliance level. Includes IEC 61000-4-3 test.

20

MANUFACTURER: Shenzhen Yimi Life Technology Co., Ltd
Add: 305, Building A, Tengbo Industrial Park, Changshangjiang Street, Longbei Village, Pingshan District, 518118, Shenzhen
Tel: +86 755-86573112
Email: hnsd@mypo2.com
Web: www.mypo2.com

EC REPRESENTATIVE
Shanghai International Holding Corp. GmbH (Europe)
Address: Eiffestrassee 80, 20537, Hamburg, Germany
E-mail: shholding@hotmail.com