

# TEST REPORT



ORIGINAL

**Test Report No.** 70.405.20.12619.01  
**Dated** 2020-05-13

Applicant : KUNSHAN JIEHONG NONWOVEN PRODUCT CO., LTD.  
Address : 895 Xinle Rd, Dianshanhu Town 215345 Kunshan  
Attn : Zhan Hao

Sample Name : Isolation gowns  
Sample Type : L/200426  
Season : -  
Manufacturer : KUNSHAN JIEHONG NONWOVEN PRODUCT CO., LTD.  
No. of samples : 13 pieces  
Supplier : -  
Buyer : -  
Country Of Origin : -  
Goods Export To : -

Care Label Provided : -

Receipt Date Of Sample : Received On 2020-05-06  
Date Of Testing : From 2020-05-06 To 2020-05-13  
Sample Submitted : The sample(s) was (were) submitted by applicant and identified.  
The product's function is declared by the applicant.  
The testing point(s) was (were) selected by the applicant.

Test Result : Refer To Next Page.

TÜV SÜD Certification and Testing (China) Co., Ltd. Shanghai Branch

Testing Center

Prepared by:

Kathy Xu  
Softlines Department

Authorized by:

Andy Wang  
Softlines Department

Note: (1) The TÜV SÜD Certification and Testing (China) Co., Ltd. "General Terms & Conditions" applied.

For full version, please visit: <http://www.tuv-sud.cn/cn-scn/terms-and-conditions>

(2) The results relate only to the Items tested, (3) The test report shall not be reproduced except in full without the written approval of the laboratory.

(4) Any use for advertising purposes must be granted in writing. This technical report may only be quoted in full. This report is the result of a single examination of the object in question and is not generally applicable evaluation of the quality of other products in regular production. For further details, please see testing and certification regulation, chapter A-3.4.

Laboratory:  
TÜV SÜD Certification and Testing  
(China) Co., Ltd. Shanghai Branch  
Testing Center Softlines Lab

Phone: +86 21 6037 6300  
Fax: +86 21 6037 6350  
E-mail: [softlines@tuv-sud.cn](mailto:softlines@tuv-sud.cn)  
[www.tuv-sud.cn/softlines](http://www.tuv-sud.cn/softlines)

Regd. Office:  
TÜV SÜD Certification and Testing (China) Co., Ltd. Shanghai  
Branch Testing Center Softlines Lab  
2F Building B, No. 1999 Du Hui Rd, Minhang District, Shanghai  
201108 P.R. China


Page 1 of 2

# TEST REPORT

Test Report No. 70.405.20.12619.01  
 Dated 2020-05-13



Description Of The Test Subject

Sample	Description	Photo
001	Isolation gowns in yellow	

## Test Results

### 1 Water Resistance-Impact Penetration Test AATCC 42:2017e

	001-front chest						
	1#	2#	3#	4#	5#	6#	
Weight of blotter gained(g)	0.0	0.0	0.0	0.0	0.0	0.0	
	001-front chest						
	7#	8#	9#	10#	11#	12#	13#
Weight of blotter gained(g)	0.0	0.0	0.0	0.0	0.0	0.0	0.0
	001-sleeve						
	1#	2#	3#	4#	5#	6#	
Weight of blotter gained(g)	0.0	0.0	0.0	0.0	0.0	0.0	
	001-sleeve						
	7#	8#	9#	10#	11#	12#	13#
Weight of blotter gained(g)	0.0	0.0	0.0	0.0	0.0	0.0	0.0

### 2 Hydrostatic head Test

AATCC 127:2017(2018)e Option 2; Hydrostatic Head; Rate Of Increase Of Water Pressure: 60cm/min., face Side Facing Water

<u>As Received</u>	001-front chest						
(cm.H <sub>2</sub> O)	1#	2#	3#	4#	5#	6#	
	133.6	109.1	120.4	116.3	123.4	127.5	
<u>As Received</u>	001-front chest						
(cm.H <sub>2</sub> O)	7#	8#	9#	10#	11#	12#	13#
	120.4	108.1	123.4	120.4	123.4	114.2	126.5
<u>As Received</u>	001-sleeve						
(cm.H <sub>2</sub> O)	1#	2#	3#	4#	5#	6#	
	130.6	99.8	115.3	118.3	119.3	121.4	
<u>As Received</u>	001-sleeve						
(cm.H <sub>2</sub> O)	7#	8#	9#	10#	11#	12#	13#
	129.5	119.3	108.1	118.3	126.5	128.5	115.3

-End of The Test Report-